

Drägerwerk AG & Co. KGaA, 23542 Lübeck

**To the customers and users with  
Carbon dioxide (CO<sub>2</sub>) / oxygen (O<sub>2</sub>)  
terminal units - DIN standard**

May 2017

**Important safety notice!!!**

**Possible gas flow on DIN CO<sub>2</sub> connectors in DIN O<sub>2</sub> terminal units**

Dear Sir/Madam,

As part of our global market monitoring we have become aware of individual cases in which temporary gas flow occurred due to improper insertion of a CO<sub>2</sub> connector into an O<sub>2</sub> terminal unit by inserting with excessive force.

Due to the temporary gas flow there is the risk that O<sub>2</sub> instead of CO<sub>2</sub> is applied to a patient during a minimally invasive procedure. There were no severe patient injuries in any of the known cases.

CO<sub>2</sub> terminal units are installed in healthcare facilities together with O<sub>2</sub> terminal units in selected functional areas such as operating theatres or endoscopy. A relevant risk of confusing the terminal units only occurs where the CO<sub>2</sub> and O<sub>2</sub> Outlets are installed in direct proximity to each other.

In case of an improper insertion of a CO<sub>2</sub> connector into an O<sub>2</sub> terminal unit, the connector does not lock and is pressed out again by the gas pressure after a short period of time.

According to the applicable standard DIN 13260-2, the O<sub>2</sub> and CO<sub>2</sub> connector geometries are so similar that a CO<sub>2</sub> connector can also be inserted into an O<sub>2</sub> terminal unit. According to DIN EN ISO 9170-1, the distinctive feature of the coloured or colour-neutral version of the coupling is only the gas type symbol. According to the instructions for use, it must be ensured during application that the marking of the connector corresponds to the marking of the release sleeve of the terminal unit so that improper insertion can be avoided.

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53-55  
23558 Lübeck, Germany  
Postal address:  
23542 Lübeck, Germany  
Tel. +49 451 882-0  
Fax +49 451 882-2080  
info@draeger.com  
www.draeger.com  
VAT no. DE135082211

Bank information:  
Commerzbank AG, Lübeck  
IBAN: DE95 2304 0022 0014 6795 00  
Swift code: COBA DE FF 230  
Sparkasse zu Lübeck  
IBAN: DE15 2305 0101 0001 0711 17  
Swift code: NOLADE21SPL

Company headquarters: Lübeck  
Commercial register:  
Municipal Court Lübeck HRB 7903 HL  
General partner: Drägerwerk Verwaltungs  
AG  
Company headquarters: Lübeck  
Commercial register:  
Municipal Court Lübeck HRB 7395 HL

Chairman of the Executive Board  
of  
Drägerwerk AG & Co. KGaA and  
Drägerwerk Verwaltungs AG:  
Prof. Dr. Nikolaus Schweickart  
Board of directors:  
Stefan Dräger (Chairman)  
Rainer Klug  
Gert-Hartwig Lescow  
Dr. Reiner Piske

Page 2/3

Based on the findings of the investigation carried out, Dräger has developed an improved component for the O<sub>2</sub> terminal unit that reliably prevents gas flow even when a current, compliant CO<sub>2</sub> connector is accidentally inserted using excessive insertion force. In addition, the instructions for use have been updated and now include a more detailed description of the insertion process.

During the review of the events reported in the field, we discovered that also CO<sub>2</sub> connectors were used that do not comply with the current standard DIN 13260-2 any more. These connectors have a smaller internal diameter of 3.5 mm which can lead to a temporary gas flow due to the tight tolerance requirements of the previous standard in connection with the plug-in coupling (even after the update). As part of an adjustment of the standard in 2004, the internal diameter of the CO<sub>2</sub> connector was changed to 5.9 mm.

Please check the internal diameters of all CO<sub>2</sub> connectors in your healthcare facility. We recommend using only CO<sub>2</sub> connectors that comply with the current standard DIN 13260-2.

Please contact Dräger in order to update the affected O<sub>2</sub> terminal units. Within the warranty period or in case of regular service by Dräger, this is done free of charge. For this purpose, please complete the attached "Customer feedback and order" form and return it to us. We will then contact you to coordinate all required activities.

Please pay attention to the following points until the upgrade is carried out:

- The gas type label of connector and terminal unit must correspond
- The coloured gas type marking (if applicable) of connector and terminal unit must correspond
- The insertion of the connector into the terminal unit must be smooth
- Pay attention to the clicking sound when the connector locks.
- Slightly pull on the connector after it has locked to check the plug connection. The connector must not loosen.
- The connector must not be pressed in when in operating position.

We regret any inconvenience caused by this and would ask you to inform all personnel in your facility accordingly. Unfortunately, this activity is unavoidable as a preventive measure to increase patient and user safety.

The competent authorities have been informed about this activity.

With many thanks in advance for your support.

Yours faithfully,

Page 3/3

Dr. V. Prapavat  
Vice President Workplace Infrastructure  
Product Management & Application Marketing  
Drägerwerk AG & Co. KGaA

Find attached:  
Customer feedback and order  
Supplementary sheet for the instructions for use