
# MANUAL RESUSCITATION BREATHING SYSTEM


en	Manual Resuscitation Breathing System
fr	Reanimateur Manuel
de	Manuelles Rettungs-Beatmungs-System
es	Sistema de Resucitación Manual
pt	Sistema Respiratorio para Ressuscitacao Manual
it	Circuiti per rianimazione autogonfianti (BVM)
nl	Manuele Resuscitatie Set
no	Resusciterings-system
fi	Käsiresuskitaattori
sv	Manuell Raddningsblasa
da	Manuelt Ventilationssystem
el	Συσκευή ανανήψης χείρος, αμπου
lt	Rankinė gaivinimo sistema
pl	System do manualnej resuscytacji
ru	Ручные реанимационные системы (мешки)

cs	Manuální resuscitační dýchací systém
hu	Kézi újraélesztő rendszer
sl	Ročni sistem za reanimacijo s predihavanjem
lv	Manuālā elpināšanas sistēma
et	Hingamiskotisüsteem manuaalseks ventilatsiooniseks
bg	Система за ръчно възстановяване на дишането
sr	Ručni sistem za reanimaciju
ro	Sistemul de resuscitare pentru ventilație manuală a plămânilor
sk	Resuscitačný dýchací system manuálny
hr	Manualni sustav za reanimaciju
tr	Elle İşletilen Solunum Canlandırma Sistemi
ja	インターサージカルBVM
zh-t	手動式復蘇呼吸系統
zh-s	手动式复苏呼吸系統
ar	نظام إنعاش تنفسي يدوي.


1.


2.





3.

4.

## TECHNICAL DATA

Resuscitation bag size	Patient weight, kg	Volume, ml	Tidal volume, ml	Pressure relief pop-off setting ( $\pm 5H_2O$ ) (* If applicable)	Maximum Dimensions (mm)	Weight (kg)	Deadspace, ml
Adult	>40	1500	610	60/40 $\pm 5H_2O$ *	322mm x 165mm	0.369	8.30ml
Small adult	10-40	1000	450	40 $\pm 5H_2O$ *	294mm x 136mm	0.295	8.30ml
Child	5-10	550	240	40 $\pm 5H_2O$	244mm x 112mm	0.210	8.30ml
Paediatric (Infant)	0-5	280	40	40 $\pm 5H_2O$	243mm x 92mm	0.187	8.30ml

**CAUTION:** Resuscitation breathing systems must only be used by persons who have received adequate training and instruction from qualified medical personnel. Training should include bag - valve - mask resuscitation on realistic training mannequins. Inappropriate use of these products may result in harm to the patient. These resuscitators cannot be prevented from entraining air from the environment.

This may be dangerous to the patient if resuscitation is being performed in a hazardous atmosphere, if in doubt use a validated filter. In the presence of high oxygen concentrations, an explosion may be caused by naked flames, or by oil and grease on the resuscitator. When used by a trained operator, Intersurgical Resuscitation Breathing Systems can be used to maintain ventilation safely and efficiently, either in situations where resuscitation is required, or where breathing is insufficient. All systems can be connected to an oxygen supply. The oxygen concentration delivered will depend upon the oxygen flow rate, the presence of a reservoir bag, and the tidal volume. Resuscitation can be performed without supplementary oxygen. Vomit can be cleared from the patient connection end by shaking and squeezing the bag.

Intersurgical Resuscitation Breathing Systems are available with 4 sizes of self-inflating resuscitation bags.

1. Check the bag, valve and oxygen delivery for function and integrity.
2. Connect the oxygen supply line, if present and clinically necessary.
3. For systems incorporating pressure control, maximum pressure achievable is indicated on device (@30L/min)

This can be manually overridden by holding it closed.

<b>O<sub>2</sub> Concentration - 12bpm (O<sub>2</sub> %)</b>			
Size	Oxygen (%)		
	5L/min	10L/min	15L/min
Adult	60	94	97
Small adult	72	94	95
Child	79	93	98
Paediatric (Infant)	92	98	98

IQR 159 ID 88682

Environmental Operating Conditions: -18°C +50°C

Storage Temperature: -40°C +60°C

Inspiratory Resistance: <3.0cm H<sub>2</sub>O at 50L/min (1cm H<sub>2</sub>O ~ 0.1kPa)

Expiratory Resistance: <3.0cm H<sub>2</sub>O at 50L/min

Resuscitator maximum weight: 0.6kg

\*Under extreme cold conditions performance may be reduced.

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Taille de ballon	Poids du patient (kg)	Volume (ml)	Volume curant inflation (ml)	Réglage de décompression d'évacuation ( $\pm 5H_2O$ ) (* si c'est approprié)	Dimensions maximal (mm)	Poids (kg)	Espace vide (ml)
Adulte	>40	1500	610	60/40 $\pm 5H_2O$ *	322mm x 165mm	0.369	8.30ml
Petit Adulte	10-40	1000	450	40 $\pm 5H_2O$ *	294mm x 136mm	0.295	8.30ml
Enfant	5-10	550	240	40 $\pm 5H_2O$	244mm x 112mm	0.210	8.30ml
Nourrisson	0-5	280	40	40 $\pm 5H_2O$	243mm x 92mm	0.187	8.30ml

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Rx ONLY 0120

**INTERSURGICAL**<sup>®</sup>  
COMPLETE RESPIRATORY SYSTEMS

Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK

T: +44 (0)118 9656 300 F: +44 (0)118 9656 356 [info@intersurgical.com](mailto:info@intersurgical.com) [www.intersurgical.com](http://www.intersurgical.com)

**Ireland**

T: +353 (0)1 697 8540  
[info@intersurgical.ie](mailto:info@intersurgical.ie)

**France**

T: +33 (0)1 48 76 72 30  
[info@intersurgical.fr](mailto:info@intersurgical.fr)

**Deutschland**

T: +49 (0)2241 311063  
[info@intersurgical.de](mailto:info@intersurgical.de)

**España**

T: +34 91 665 73 15  
[info@intersurgical-es.com](mailto:info@intersurgical-es.com)

**Portugal**

T: +351 219 108 550  
[info@intersurgical.pt](mailto:info@intersurgical.pt)

**Italia**

T: +39 0535 20836  
[info@intersurgical.it](mailto:info@intersurgical.it)

**Nederland**

T: +31 (0)413 243 860  
[info@intersurgical.nl](mailto:info@intersurgical.nl)

**Sverige**

T: +46 08 514 30 600  
[info@intersurgical.se](mailto:info@intersurgical.se)

**Danmark**

T: +45 28 69 80 00  
[info@intersurgical.dk](mailto:info@intersurgical.dk)

**Lietuva**

T: +370 387 66611  
[info@intersurgical.lt](mailto:info@intersurgical.lt)

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T: +7 495 7716809  
[info@intersurgical.ru](mailto:info@intersurgical.ru)

**Česká Republika**

T: +420 272 911 814  
[info@intersurgical.cz](mailto:info@intersurgical.cz)

**Türkiye**

T: +90 216 468 88 28  
[info@intersurgical.com.tr](mailto:info@intersurgical.com.tr)

**South Africa**

T: +27 (0)11 444 7968  
[info@intersurgical.co.za](mailto:info@intersurgical.co.za)

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T: +86 519 69817000  
[info@intersurgical-cn.com](mailto:info@intersurgical-cn.com)

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T: +81 (0)3 6809 2046  
[info@intersurgical.co.jp](mailto:info@intersurgical.co.jp)

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T: +886 4 2380 5430  
[info@intersurgical.com.tw](mailto:info@intersurgical.com.tw)

**Pilipinas**

T: +63 2820 4124  
[info@intersurgical.ph](mailto:info@intersurgical.ph)

**USA**

T: +1 800 828 9633  
[support@intersurgicalinc.com](mailto:support@intersurgicalinc.com)

**Colombia**

T: +57 1 742 5161  
[info@intersurgical.co](mailto:info@intersurgical.co)

**Australia**

T: +61 (0)2 8048 3300  
[info@intersurgical.com.au](mailto:info@intersurgical.com.au)

Distributed in the USA by Intersurgical Incorporated, 6757 Kinne Street, East Syracuse, NY 13057

T: (800) 828-9633 F: (315) 451-3696 [support@intersurgicalinc.com](mailto:support@intersurgicalinc.com) [www.intersurgical.com](http://www.intersurgical.com)

7762 Issue 8 IP 08.16