



Urgent Field Safety Notice (Removal)

Cordis Vista Brite Tip® & ADROIT™ Guiding Catheters

February 07, 2019

Dear Valued Customer,

The purpose of this communication is to inform you Cordis is recalling 173 lots of Cordis Vista Brite Tip® and ADROIT™ Guiding Catheter product.

| | |
|---|--|
| <p>Recall Overview:</p> | <p>Cordis has initiated a recall for 173 lots of Vista Brite Tip® and ADROIT™ Guiding Catheter due to frayed pieces of the mounting card being inside the primary packaging. See below image.</p>  <p>Frayed pieces of the mounting card were found in the sterile primary package. While the material may be identified when the product is opened and prepared for use, there is a potential for harm if the material is not discovered or if the device is not prepared properly. Prior to and during preparation of the device, if the frayed material is identified, the user would be prompted to exchange the device for another one, resulting in a pre-procedural or intra-procedural delay. However, if undiscovered, the frayed material could enter the patient's vasculature potentially resulting in ischemia, necrosis, or the need for additional intervention. The sterility of the product is not affected. Cordis has not received any reports of patient harm related to this issue.</p> <p>Usage</p> <p>The Vista Brite Tip® and ADROIT™ Guiding Catheter is intended for use for intravascular introduction of interventional/diagnostic devices into the coronary or peripheral vascular systems.</p> |
| <p>Details on Affected Devices, to assist in identification of the product involved:</p> | <p>Identification</p> <p>See Table 1 below for catalog number and lot number list</p> <p>Example labeling:</p>  |

| | |
|--|--|
| Why you are being contacted: | You are receiving this letter because our records indicate that you have purchased impacted Vista Brite Tip® and ADROIT™ Guiding Catheter lot numbers. |
| Actions requested on your part: | <ol style="list-style-type: none"> 1) Read this Field Safety Notice (Removal) letter. 2) Immediately check your inventory to confirm whether or not you have any units from the affected lots in your possession. Identify and set aside any units from the affected lots in a manner that ensures the affected product will not be used. Check all storage and usage locations. 3) Review, complete, sign and return the enclosed Acknowledgement Form in accordance with the directions on the form. 4) Return all affected product to the Cardinal Health distribution center. Please contact your local sales representative to facilitate return of the affected product, if necessary. Your sales representative will inform you of the product replacement or credit options. 5) Share this letter with others in your facility who need to be made aware of this recall and with any other facility that may have been sent the affected units of Vista Brite Tip® and ADROIT™ guiding catheter product from your facility. If any units of the affected lots are found to be at the other facility, please arrange the return of the units. 6) Maintain awareness of this notice until all affected product has been returned to Cordis and keep a copy of this notice with the affected product. |
| Description of the problem: | <p><u>What is the issue?</u> Cordis became aware of frayed pieces of the mounting card being inside the device primary packaging.</p> <p><u>Why are we recalling this product?</u> Frayed pieces of the mounting card were found in the sterile primary package. While the material may be identified when the product is opened and prepared for use, there is a potential for harm if the material is not discovered or if the device is not prepared properly. Prior to and during preparation of the device, if the frayed material is identified, the user would be prompted to exchange the device for another one, resulting in a pre-procedural or intra-procedural delay. However, if undiscovered, the frayed material could enter the patient's vasculature potentially resulting in ischemia, necrosis, or the need for additional intervention. The sterility of the product is not affected.</p> <p><u>What other actions is Cordis taking?</u> Cordis is investigating the root cause and will take appropriate corrective action. Cordis has not identified any other lots that may be affected.</p> |
| Available Assistance: | If you have any questions regarding this recall, please contact your local sales representative or local sales office. |
| Additional Information: | <p><u>Regulatory Notification</u> The applicable regulatory agencies and notified body are being notified that Cordis is voluntarily taking this action.</p> |

We apologize for any inconvenience this communication may cause. We know that you place high value in our products, and we appreciate your cooperation in this matter. Cordis is committed to maintaining your confidence in the safety and quality of the products that Cordis supplies.

Respectfully yours,

A handwritten signature in black ink, appearing to read "Miguel Ávila". The signature is fluid and cursive, with a long horizontal stroke at the end.

Miguel Ávila
Vice President, Global Quality and Regulatory Affairs
Cordis Corporation

Table 1

| Catalog # | Lot Number | Catalog Number | Lot Number |
|-----------|--|----------------|--|
| 67019055 | 17744955 17750783 17757121 | 588840P | 17756076 17816499 |
| 67021055 | 17735533 17735534 17744429 17744430 17750785 17757116 | 598943P | 17746147 17807995 17811028 17814511 17818155 |
| 67021255 | 17748743 17750476 17755096 17757122 | 588843P | 17777712 17816501 |
| 67021455 | 17753002 | 588845P | 17729208 |
| 67027055 | 17779364 | 588841P | 17753948 |
| 67027855 | 17815946 | 588857P | 17756074 17816500 |
| 67028055 | 17805148 17815949 | 588858P | 17756075 |
| 77821055 | 17746022 17762925 | 598945P | 17755015 |
| 77821255 | 17733016 17747583 17749800 17752698 17767925 17769354 17771053 17772929 | G780GOND | 17724819 17742581 17756975 17771328 17794750 17814264 17816167 |
| 77822455 | 17722360 17735189 | SM7673 | 17750494 17817649 |
| 77827055 | 17733019 17746026 17747584 17749793 17752703 17754875 17756762 17762928 17771058 17772931 17816029 17817637 | 588846P | 17733005 17736643 17737608 17742559 17740605 17743099 17743101 17745965 |
| 77827855 | 17722359 17740611 17778655 | 77828055 | 17724480 17767923 17771049 17772926 |

| | | | |
|----------|----------|----------|----------|
| 67221255 | 17753959 | 6720540E | 17800760 |
|----------|----------|----------|----------|

Table 1 (continued)

| Catalog # | Lot Number | |
|-----------|------------|-------------------------------|
| 55805400 | 17801566 | 17806219 17808004 17809486 |
| 67005400 | 17802170 | 17811837 |
| | 17802171 | 17811839 |
| | 17803141 | 17813158 |
| | 17803142 | 17813159 |
| | 17803143 | 17815025 |
| | 17803536 | 17817001 |
| | 17803539 | 17817002 |
| | 17803540 | 17817004 |
| | 17804225 | 17817285 |
| | 17804226 | 17818563 |
| | 17804764 | 17818564 |
| | 17805447 | 17819221 |
| | 17805899 | 17819222 |
| | 17805901 | 17821534 |
| | 17806330 | 17800613 |
| | 17806331 | 17800614 |
| | 17806332 | 17800615 |
| | 17807603 | 17802172 |
| | 17807626 | 17803537 |
| | 17807628 | 17803538 |
| | 17807629 | 17804026 |
| | 17809536 | 17804222 |
| | 17809537 | 17804223 |
| | 17809538 | 17804224 |
| | 17809786 | 17804765 |
| | 17809791 | 17804766 |
| | 17809792 | 17805898 |
| | 17809793 | 17805900 |
| | 17810429 | 17805902 |
| | 17810430 | 17807627 |
| | 17810432 | 17811835 |
| | 17810433 | 17815021 |
| | 17811219 | 17815022 |
| | 17811220 | 17815023 |
| | 17811221 | 17815024 |
| | 17811834 | 17817003 |
| | 17811836 | 17819220 |
| 67205400 | 17806706 | |
| | 17807595 | |
| | 17811677 | |
| | 17812321 | |
| | 17814278 | |
| | 17814279 | |
| | 17814640 | |
| | 17814641 | |
| | 17814642 | |
| | 17814643 | |
| | 17815340 | |
| | 17806214 | |