

Urgent Field Safety Notice

| | |
|--------------------------------|---|
| Affected Product | ONLY CERTAIN CODES - Note: see Addendum "A" for full details. Feeding Tubes, Nelaton Catheters, Levins Tubes, Rectal Catheter, Ryles Tubes, Suction Catheters, Stomach Tubes, Medical Packs. |
| Type of action | Voluntary Recall / Product Disposal |
| Pennine Healthcare Ref: | PHFSN0219 |
| Product Code(s) | Various – As per Addendum "A" |
| Lot Number | Various – As per Addendum "A" An example in this range of LOT numbers is: "06D15" |

8th May 2019

Dear Customer,

Description of the problem:

This communication is to inform you that Pennine Healthcare has voluntarily initiated a recall to remove multiple LOTs of the products detailed within this Field Safety Notice (FSN).

Through post production activities we have identified that there is a small possibility that the products detailed in this FSN may not have been manufactured with the correct eyes (holes) punched into the tube. The defect in the products should be visible to a trained user, however, if the user missed the defect, it would mean the products would not perform to the required expectations. There is a potential for serious harm should an alternative non-defective device not be available at the time.

Actions to be taken by the distributor / user:

1. Identify and quarantine all stock of the affected batch / LOT numbers immediately.
2. If you are a distributor you are required to confirm to Pennine that you have provided this FSN to all your customers and contained any returned stock.
3. Complete the attached Reply Form (Addendum B) to confirm that you have read and understood the contents of this Field Safety Notice.
4. Upon receipt of the completed reply form, Pennine Healthcare will arrange for the return or confirm how destruction of the affected product should be organised. (Proof of destruction and records of type and quantities involved will be needed).
5. For further clarification of quarantined or affected stock replacements or credits, please contact your local distributor or if supplied directly from Pennine our customer services using:

- UK markets – customer.services@penninehealthcare.co.uk tel. +44 (0)1332 794 880
 - International Markets - exports@penninehealthcare.co.uk tel. +44 (0)1332 794 880
6. Maintain awareness of this Field Safety Notice until all stock of the affected LOT batch has been destroyed or returned to Pennine Healthcare.

Identification of affected products:

Use the Product Codes and LOT numbers as detailed in Addendum A to identify the affected devices.

Transmission of this Field Safety Notice:

This notice should be passed on to all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred.

Please maintain awareness on this notice until all required actions have been performed within your organisation.

We confirm that this field safety corrective action has been communicated to the relevant competent authorities and distributors.

Pennine is committed to providing quality products to our customers and we sincerely apologise for any inconvenience this notice may cause.

For and on behalf of Pennine Healthcare:



Simon Motlib
Quality Assurance Manager

Addendum A - Affected Medical Devices – (NB - Only these product codes are affected)

| Product Code | Batch Numbers |
|---------------------------|--|
| Medical Pack FRP-04400 | 02C16,09C16,10M15,15D16,17B16, |
| Medical Pack FRP-04700 | 05K18,06L18,12B18,19M18,29K18, |
| Medical Pack FRP-04701 | 04B19,14C19,17D19,23A19, |
| Feeding Tubes (FT) | |
| FT-1608/105 | 01C17,01M17,06C18,07L16,09E18,09L17,11K18,12J17,18J14, |
| FT-1610/105 | 09F14,11H17,15J15,17G15,17L14,28A19,30K17, |
| FT-1908/105 | 03J14,04H16,04J15,06H15,08A19,09B15,09F14,09G15,10G17,11B19,11G14,11L15,13E15,15D16,16K17,16L18,17D18,18C19,18K18,21L17,22M14,25G17,27E14,27L18,29H17,30G14, |
| FT-1908/60 | 01C16,02G15,04A18,04D19,04M18,06K14,12C18,13L18,20F17,20J18,25K18,26D16, |
| FT-1910/105 | 03G14,18J14, |
| FT-1910/60 | 04M18,10A19,20J18,25K18, |
| Levins Tubes (LT) | |
| LT-2208 | 01C17,01E15,01M17,02B18,04G18,05F15,07H14,07L16,09F17,10G17,11H17,13F16,16B15,17F15,21A15,22H16,22K18,25E18,28C18,28K14,30A15,30D18,31C15, |
| LT-2210 | 01B17,01E15,02M16,04G17,05F15,09F17,11H17,12K16,13F18,14D15,15A15,17D18,17F15,17H18,18B16,19C15,23E16,24L17,26E15,26F14,27B18,27C17,28A19,29A16,29G14,30D18, |
| LT-2212 | 01B17,01E15,01M17,03E17,03K17,04B16,05F15,06D18,09F14,12M14,15A15,15A18,16B15,16F15,16M16,17B16,17K14,18L14,18M18,19H15,19M16,22A16,22D16,22E15,22F17,22H16,23D18,23E16,25E18,26K15,27G18,29B16,30D18,31C17,18J14, |
| LT-2214 | 01B17,01C17,02D19,02M16,03G14,03K16,04B16,04E18,06H15,06L14,06M17,07F18,08F15,09E14,10D17,10L16,11B16,12E15,12M14,13A16,13K17,19E14,19F18,20K14,20K15,21A15,21D16,22H14,22M16,24A18,25F14,25J18,26E15,26J17,26K15,27K16,27L14,29F17,29K18,31E16, |
| LT-2216 | 01J16,01L18,02D19,02E14,03E17,03G17,03H17,03K17,04J15,05A16,05B19,05M16,06C17,06C19,06D18,06K15,06M17,07L16,11B16,11D18,11G14,11J14,11K18,12A18,12E14,13F16,14A19,14C18,14D15,15A15,15D16,16F17,16L16,16L18,17B16,17E16,17E18,17L14,19C15,19E14,19J14,19J16,19J17,22C19,22F15,22G14,22G16,22M14,23G15,24B16,26E15,27E14,27F18,27G17,27G18,28B17,28C17,28K15,29B16,30A18,30D18,30K17, |
| LT-2218 | 01K18,01M17,02E16,02M16,03K17,04E18,05A16,06C17,06D17,09G15,10B16,10C15,10D19,10G17,11D16,12E14,13F18,13K15,14C18,15A15,15B17,15J17,16F15,16G18,16L18,17E16,17L16,18H14,18L15,19J14,20D15,21A19,22F17,22G14,22M14,22M16,24A18,26E15,26H16,27L14,28H15,28K15,29F16,29K18,30G18, |

| | |
|-------------------------------|--|
| LT-2220 | 01J16,02G14,03E18,06A17,08F18,09K17,13G16,17A17,17G17,19B15,27F15,29E15,31C17, |
| LT-2408 | 04B16,04H16,04J15,05D16,08A19,11G14,11H16,12J17,18H16,18J14,18L14,22G15,25F14,25J14,26J17,30J16, |
| LT-2410 | 03G14,06C18,08A19,08D16,09F14,15K18,16J16,17G15,17L14,18J14,23G15,29A16,30K17, |
| LT-2412 | 01M17,03G14,04J15,06C18,06E14,08D16,11H16,14D15,15K18,16G15,16M16,17L14,18G18,19J14,20C18,24J15,26J16,27C17,29B16,30D18,31C17, |
| LT-2414 | 02E14,03D17,03G14,03M15,04B16,04H16,05D16,05K15,08J16,09F17,11H16,14C18,14D15,14M17,15C17,15D16,15J15,15K18,16F14,16G15,17H17,17L14,19G18,19J14,22M16,25G17,25J14,27K16,30A15,30E17, |
| LT-2416 | 02E14,03G14,04H16,05D16,06C18,06D17,06K15,07L16,08M15,10B15,11A17,11H16,11K18,14A19,14D15,15D16,16F14,16F17,16J15,17H17,18L14,19J14,19J16,20M17,22E17,22G15,22M14,23G18,24B16,26J16,27C17,27G17,29F16,30D18,30K17, |
| LT-2418 | 03A18,03H17,03M15,04B16,05F14,08A19,09B15,10C15,10G17,11D16,11H17,12A18,14D15,15D16,15G14,15K18,16J15,17D18,17L14,19H15,19J14,20K16,22G15,22M16,25G17,26C18,28G16,29B16,30E17,30J16,31C17, |
| LT-2420 | 02G14,02L17,03B15,10L14,16J16,18J14, |
| Nelaton Catheters (NC) | |
| NC-1220/SW | 02J16,03F16,03G14,03K14,04A19,04D17,04G18,06B17,07C18,07D15,07L17,08E18,08F18,08L16,09A15,09L15,10D18,10G15,11K17,16D15,16G18,16L15,17E17,18B19,18K18,20J17,20L14,21J16,22H16,22H18,22J14,24B15,24G14,25E16,26D18,28J18,29E15,30K15, |
| NC-1220/SW/1 | 04A17,06B17,07D15,10C17,11J18,20L14,24B16,25J15,26K15, |
| Rectal Catheters (RC) | |
| RC-1820 | 02C16,02J16,03F16,04A17,04A19,04D17,05F14,06B17,07C18,07D15,08F18,10D18,10M18,12D16,13M17,15D19,16G18,16M14,17E17,18B19,19K17,20G17,22H18,23B17,24B15,24F14,25H15,26K15,26K16,28J18,29F15, |
| RC-1820/1 | 01H14,06E15,07C18,07K15,09A15,09H17,10D18,11A18,11J18,12D16,16L15,17E17,18B19,20G17,22H16,22H18,23B17,26K16,28H14,28J18,29F15,29F16,30L17, |
| Ryles Tubes (RT) | |
| RT-2008 | 01M17,02B18,03M14,04B16,04H16,04J15,05C15,05E15,08A15,08B17,08J16,09E18,09L17,11B19,11H17,12E14,13E15,14A19,15B17,16B15,16L18,17H18,18C19,18G18,18H16,19G17,20D15,21L17,22B18,22K18,23E16,26H16,26J17,27E14,28K14,28K15,31C15,31E16, |
| RT-2008/BULK | 06L15,16C17,21C16, |
| RT-2008/L | 04G18,09F14,11K18,17K14,22H16,22K18,28C18, |
| RT-2010 | 01E15,02B18,02M16,03J14,04G17,04J15,05A16,05F15,06C18,06K15,06L14,07B17,07L16,08A15,09C15,09E16,09F14,09F17,10D17,11A17,11H17,13B18,13F18,13M18,14D15,14E15,14J15,15K18,16J16,16L16,16L18,17D18,17G15,17H18,18C16,18C19,18L15,19E14,19G17,20M17,21F16,23B15,24J18,24L17,25C15,25E18,25G17,26H16,27C17,27K16,27L14,27L18,28A19,28K14,29G14,30A15,30G18,31E16, |
| RT-2010/BULK | 10A17,15C16,25J17,27C19,14M15, |

| | |
|--------------|--|
| RT-2010/L | 03L15,04G17,04H15,05F15,06K15,08D16,08K18,09B15,09E16,09F17,13B18,13E15,13F16,13F18,13M18,15D15,15J17,16K14,16L16,17L14,18H16,19G17,19H15,20B17,22F15,22K18,24A18,25E18,26C18,26H16,29A16,31G14, |
| RT-2012 | 01L18,01M17,03E17,03F14,04B16,04H15,05C15,06C17,06C18,06D17,06E14,06K15,06L14,07C16,07H14,08A15,08J16,08M15,09E16,09E17,09F14,09F17,10D17,10L16,11D18,11G18,11K17,11K18,11L15,12K16,13E15,13F16,13F18,13G16,13J18,14A19,14K15,15A18,15J17,16G15,16K14,16L16,16M16,17A17,17H17,18D19,18G18,18H16,18L14,18M18,19G17,19H15,20B17,20C18,20D15,20M17,21L17,22B18,22C19,22D16,22F17,22H14,22J14,22K18,22L18,23D18,23K17,24A18,24H18,24J18,24L16,24L17,25E18,26C19,26J16,26K15,27C17,27E15,27G17,27H15,27L14,28A19,29C16,29F15,29G14,30G18,31C15,02K14, |
| RT-2012/BULK | 02C17,06L15,12A18,21C16, |
| RT-2012/L | 05J17,07C16,08A15,08K18,09E16,09F14,11D16,11G18,11H16,13E15,13F16,13F18,13J18,14A19,15D15,15J17,16F15,16G15,16L16,18B16,18D19,18M18,19H15,20B17,21F16,21L17,22L18,23E16,24H18,24J15,24J18,27E14,27G17,30D18, |
| RT-2014 | 01D15,01J16,01K18,02B18,02E14,02F14,03A18,03H17,03J18,03K16,04D19,04E18,04G18,04J14,05H15,05M16,06C15,06M17,07F18,07H14,07K15,08A15,08A19,08F15,09G15,09L18,11B19,11D17,12E14,12E17,12F17,12L15,12M14,13B18,13E15,13K17,13L17,14C18,14F16,14K16,14M17,15A15,15B17,15C17,15J15,15J17,15M15,16F14,16K17,16L16,16L18,16M16,17A17,17E18,17H17,18C16,18H14,19G18,19H15,19J14,20D15,20F18,20K14,20K15,21A19,21D16,22K18,23K17,24A18,25F14,25G17,25J18,26H16,26J17,26L15,27B18,27K16,27L18,28K15,29A16,29B16,29F15,29F17,29G14,29K18,29L18,30A15,31C17,31K17, |
| RT-2014/BULK | 03J15,08F17,10D17,11C16,14J17,15F15,17H16,23D15,24B16,27E16, |
| RT-2014/L | 01K18,02D19,02F14,03G14,03H17,03K16,05A16,05D16,05K15,07L18,09C15,09E18,11H16,11J17,13L17,14M17,15A15,15B17,15J17,16M16,17A17,17H18,17L14,19J14,20D15,22K18,25H15,26K15,27B18,27K16,27L14,29A16,29B16,29F17,29L18, |
| RT-2016 | 01E14,01E15,01J16,01L18,02F14,03E17,03G14,03H15,03H17,03J18,03K17,03L15,03M14,04G17,05B19,05C15,05M16,06C17,06C18,06C19,06D17,06D18,07H14,07L16,08A15,08A19,08F15,08M15,09G15,09M14,10B16,10C15,11D16,11D18,11G18,11H16,11J17,12A18,12E14,13M18,14B18,16F17,16J15,16K17,16M16,17E18,17F15,18C19,18H14,18H16,18K18,18L14,19H15,19J17,20C18,20D15,20K14,20L18,20M17,21L17,22E15,22G14,22M14,23E17,23G18,24J15,24L16,24L17,25C15,25F14,25H15,25J18,26A17,27F18,27G17,27G18,28B17,28C17,28K14,29F16,30A15,30A18,30K17,31E16,03J14, |
| RT-2016/L | 01E15,02D19,02J14,03J18,04G17,04H16,08A19,09E16,10B15,11G18,11K18,11L15,12A18,12K16,13A16,13F16,13M18,14A19,14B18,14H14,15A15,16K17,17E18,17K14,17L14,19H15,19J17,20K15,21L17,22F17,23G15,24L16,25J18,26A17,27F18,28B17,29C16,15M18, |
| RT-2018 | 01C17,01D15,01E15,01F18,01M17,02B18,02M16,03H17,03J18,03L15,03M14,04B16,07B17,07C16,08A19,08E14,08F15,09L17,10B16,10G17,10L16,11B19,11D17,11H17,12A18,12E17,12F14,13F16,13F18,13G18,13K17,13M18,14C18,15M15,16L18,17C16,17D18,17E16,17H18,17L16,18H14,18K18,18L15,21A15,21A19,21D16,21G15,22E15,22H16,22M14,24A18,25C15,25F14,25G17,26C19,26H16,27B18,27H15,27L14,28C18,28K14,29K18,30G18,31C17,02K14,11J14, |

| | |
|------------|--|
| RT-2018/L | 01E14,03H15,03M14,07B17,08K18,09L17,13E15,17H18,17L16,19H15,20K15,20K16,24A18,26C19,26J17,26K15,29B16, |
| RT-2020 | 02B17,03B15,03C16,03F14,05K15,06M18,07C19,09K17,10D18,10J18,12B16,13G18,13K14,14D15,15A16,16A17,16K17,17H15,19J17,20F14,26K16,26L15,27F15,30A18,30D15,30H16,31C15, |
| RT-2108 | 01M17,02B18,02D19,02E16,03H15,04G18,04J15,06E14,07C16,08A19,09B15,09C15,09E18,09G15,09L17,10L16,11A17,11B19,11K18,11L15,12F17,12J17,13F16,14A19,15D16,17D18,18G18,18L14,19E14,19G17,20K15,21A15,25E18,25G17,26A17,26J17,27K16,27L18,28C18,29F16,30D18,31E16, |
| RT-2108/EF | 02D19, |
| RT-2110 | 02B18,04G17,05A16,05F15,06C18,06L14,07C16,08A15,09B15,09E17,09F14,09F17,09L17,11A17,12A18,12E15,13B18,13F18,13G18,15D15,15D16,15J17,15M15,16F14,16L18,17D18,17F15,18B16,18C19,19E14,19G17,19H15,20K15,20K16,20M17,21L17,22G15,22M14,24L17,25J14,26C15,26C18,26E15,26J17,26L15,27C17,27G17,27K16,27L14,27L18,28A19,28B17,30D18,31E16, |
| RT-2110/EF | 28C19, |
| RT-2112 | 03H15,03M15,04J15,05J17,06C15,06D17,06D18,06E14,06K15,06L14,06M17,07F16,07H14,07L18,08A15,09B15,09E16,09E17,09F14,09F17,10L16,11A17,11D16,11G18,11J14,11K17,11L15,12E15,13A16,13G16,13J18,13M18,14K15,15A18,15C17,15K18,16F14,16G15,16J16,16L16,16M15,17A17,18G18,18H14,19F18,19G17,19H15,20B17,20C18,20K16,20M17,21A19,21L17,22A16,22B18,22C19,22D16,22F17,22H16,22L18,22M14,23D18,24A18,24H18,24J15,24J18,25E18,26J16,27C17,27E14,27E15,27H15,28A19,28B17,28L14,29B16,29F15,29G14,29H17,30D15,30D18,31K17, |
| RT-2114 | 01B17,02B18,02C17,02E16,03A18,03D17,03J14,03J18,03K17,03M14,03M15,04D19,04E18,04H16,04J15,05B19,05H15,05J17,05K15,06E14,06M17,07C16,07F18,07L16,08F15,09C15,09F17,09L17,09M15,10B16,11D17,11G18,11J17,11L15,12E15,12E17,12M14,13A16,13B18,13F16,13K17,14A19,14D15,14J15,14K16,14M17,15C17,15M14,16F15,16F17,16G15,16J16,16K14,16L18,17E18,17K14,18C16,18C19,18H14,19G18,19J17,20F18,20K15,21F16,22F15,22K18,22M14,23D18,23K17,24A18,24L16,25F14,25G17,25H15,25J14,25J18,26C15,26J17,27B18,27E14,27G18,27L18,28C18,28L14,29F15,29F16,29F17,29G14,29H17,29L18,30A15,31K17, |
| RT-2116 | 01C16,01F18,02C17,02F14,02K14,03E17,03J18,03K17,03L15,03M14,04G17,05A16,05B19,06C17,06C19,06D17,06D18,07F16,07H14,07M17,09C15,09M14,10B15,10G15,10L16,11A17,11G18,11J17,11K18,12M14,14A16,14A19,14B18,14J15,15B17,15D15,16F17,16K17,16L18,17D18,17E18,17F15,18H16,19J17,20C14,20C18,20K16,20L18,21F16,21L17,22C19,22F15,22H16,23E17,24A18,25H15,25J18,26C15,26K15,27F18,27G18,27H15,28B17,28L14,29H17,29L18,30A18,30D18, |
| RT-2118 | 01B17,01C17,01E14,01K18,01M17,02B18,02E16,03A18,03F14,03G14,03J18,03M15,05F15,05H15,07F16,08K18,09L17,11B19,11K17,12A18,12M14,13A16,13E14,13K15,13M18,14A19,14C18,14J15,15G14,15J17,16B15,16G18,16K14,16L18,17C16,17K14,19G17,21A15,21A19,22B18,22F17,22H16,23K17,24A18,26C15,27K16,27L18,28B17,28C18,28C19,28L14,29F16, |
| RT-2308 | 01D15,02E16,06K15,08J16,09G15,09L17,11K18,18H16,18L15,19G17,23E16,26H16,29B16,29H17,30A18,31E16, |

| | |
|-------------------------------|--|
| RT-2308/EF | 08A19, |
| RT-2310 | 02E16,03J14,04G17,04H15,06C18,07H14,08A19,08M15,09E18,09L17,09L18,11A17,11H17,13F16,13F18,15D15,16B15,17D18,17H18,19G17,19J16,22E17,22G16,22H14,23E16,24A18,27C17,27G17,28C19,30D15, |
| RT-2310/EF | 28C19, |
| RT-2312 | 01B17,03F14,04H15,05J17,07L18,09F17,10G18,11D18,11H16,11K17,13F16,13F18,15A18,15J17,15K18,16L16,18D19,18G18,20M17,22B18,22D16,22F17,22H14,22K18,22L18,23E16,24A18,24H18,25E18,25F14,26C19,27C17,28A19,29H17, |
| RT-2314 | 01C17,02B18,04D19,04H16,04J14,05A16,05H15,05M16,07F18,08F15,09E16,09L17,09L18,10B16,11G18,14K16,14M17,15B17,15C17,16F17,16G15,16L16,16M16,17A17,17H18,18H16,19G18,19J17,20C18,20D15,20F18,20K14,20K15,21A15,21A19,21F16,22H14,23B15,23D18,24A18,25G17,26L15,27L14,28C18,28G16,29C16,29F17,29H17,30E17,31C15,31C17, |
| RT-2314/EF | 20F18,27G18,28C18, |
| RT-2316 | 01B17,01J16,03H15,03J14,03J18,03K17,04G17,05A16,08A15,09L18,10B15,11D16,11G18,11H16,11J17,13M18,14H18,17B16,17E16,17E18,18K18,19C15,22B18,22E17,22G16,24L16,27C17,27G17,28A19,28B17,30D18, |
| RT-2316/EF | 02D19,06D18,12K17,23G18, |
| RT-2318 | 03H15,03K17,10G17,11D17,11L15,12E17,12K16,13K17,14H18,15B17,15D16,16G18,16L18,17D18,17E16,21A19,22F17,28G16,29H17,30E17,31C17, |
| RT-2318/EF | 04G18,26C18,27G18, |
| RT-2320 | 19J17,22B17,29H17,29K18, |
| Suction Catheters (SC) | |
| SC-1020 | 01E14,01H14,02C16,03B17,04G18,05F17,06C15,06E16,07C18,07K15,07L14,08E18,09A17,09L15,10D18,10M18,11A18,11H17,11K17,12B15,12D16,12D19,12K16,13D17,14C17,15J16,16L18,17H17,18K18,19G17,20F16,24B16,24G14,25A16,25F14,25H15,26D18,26F17,27L14,28J18,29F15,29L16,30K18,30L15,31G15,31H17,14K14, |
| Stomach Tubes (ST) | |
| ST-2408 | 01E15,02E16,04G18,09M15,10D17,11H17,12F17,16B15,17D18,22K18, |
| ST-2410 | 02B18,02K14,02M16,04F14,04G17,05C15,07H14,08D16,09E18,09F17,11H17,12J17,15M15,16G15,17C16,17H18,18B16,19G17,20K15,20M17,22E17,25G17,28A19,28C19,30G18,30K17, |
| ST-2412 | 03H15,06C17,06D17,10C15,12M14,14L14,20B17,22A16,23D18,24A18,24L16,25F14,26K15,28A19,31K17, |
| ST-2414 | 04J14,05M16,07K15,09E14,10D17,11D16,14C18,15C17,15J15,16M16,19F14,19F18,20D15,21A15,23D18,29A19,31E16,31K17, |
| ST-2416 | 06C18,06D17,10B15,11A17,11D16,11G18,11K18,14A19,14D15,16L16,20K14,22C19,22J14,22M14,26C15,29A19,29F16,30K17, |
| ST-2418 | 03H17,03M15,04G18,06K15,07F18,09A19,09B15,09F14,10B15,13M18,14C18,18K18,20D15,21A15,21A19,22F17,22H16,22M16,26C19,26J17,27B18,27E14,27G18,29F16,30K17,31C17, |

| | |
|---------------------------------|--|
| ST-2420 | 02M16,03A19,03D17,03G18,05C15,05K15,06C18,09K17,10L14,11B15,14D15,14F16,18M15,23F17,25A18,31C15, |
| X-ray Levins Tubes (XLT) | |
| XLT-2208 | 01C17,01E15,02D19,03H15,04G18,05A16,05D16,06L14,07C16,09E18,09M15,10D17,10G17,10L16,11A17,11B19,11G14,11L15,12E14,12J17,14F16,16B15,16K14,18G18,19E14,21L17,22B18,23E16,24L16,25E18,25F14,25G17,25H15,26A17,28C18,28K14,29H17,30G14,31E16, |
| XLT-2210 | 02B18,03G14,04G17,05A16,06C18,06D17,07C16,08D16,08K18,09A19,09B15,09C15,09F14,09L18,11A17,12E14,12E17,13B18,13G16,13G18,13M18,14E15,14L14,15K18,16B15,16K14,16L16,17F15,18C19,18J14,18L15,19E14,19F18,19G17,20B17,20K15,22D16,22G16,22K18,23E16,23K17,24A18,24L17,25E18,25H15,26C18,26E15,26H16,26J17,27G17,27K16,27L18,28B17,28K14,29A16,29G14,29H17,30D14, |
| XLT-2212 | 01B17,03K17,04J15,05F15,05J17,06D18,06E14,06L14,07C16,07H14,07L18,08K18,09B15,09M15,10D17,11D16,12J18,12M14,13A16,13F18,13G16,14D15,14E15,15A18,16B15,16G15,16L16,16M15,16M16,17K14,18G18,18J14,18M18,19G17,19H15,20B17,20C18,20K16,20M17,22B18,22C19,22F17,22H16,22K18,23E16,23K17,24H18,24J18,25E18,26K15,27E14,27H15,29G14,31C17,31E16, |
| XLT-2212/BULK | 11L14,11M14,19E14, |
| XLT-2214 | 01C17,01F18,02B18,02D19,02E14,02M16,03A18,03H15,03J18,03K16,04E18,05B19,05J17,06H15,06M17,07F18,07K15,08A15,08A19,09E17,09E18,09F14,10B16,10D17,11G18,11H16,12F17,12L15,13A16,13K17,14K16,14M17,15B17,15C17,15D16,15J15,15K18,15M15,16G15,16L18,17A17,17E16,17L14,18C19,18H14,18L15,18M18,19E14,19F18,19J17,20C18,20K14,21D16,21F16,22G16,23B15,23K17,24A18,24J15,25G17,25J14,26E15,26H16,26J17,26K15,26L15,27G18,27L14,27L18,29A16,29B16,29F16,29F17,29G14,29L18,30A15,30D18,31C15, |
| XLT-2216 | 01B17,01F18,01J16,01L18,02F14,03E17,03G14,03J14,03K17,04G17,04H15,04J15,05B19,05M16,06C17,06C19,07C16,07F16,07H14,07M17,08A19,08F15,09C15,09G15,10B16,11D16,11G14,11G18,11H16,11J14,11J17,11K18,11L15,12A18,12K16,12M14,13A16,13G16,13K15,13M18,14B18,14D15,14E15,14H18,15A15,15M15,16F14,16F17,16J15,16K17,16L16,16L18,17E16,17E18,18K18,18L14,19E14,19J17,20C18,20D15,20K14,20M17,21F16,22B18,22C19,22F15,22G14,22M14,23B15,23G18,24L16,24L17,25J18,26A17,27F18,27G17,27L14,28C17,28K15,29H17,30A15,30A18,30D18,30E17,30K17,31E16, |
| XLT-2218 | 01B17,01K18,01M17,02K14,02M16,03A18,03J18,04E18,05A16,05B19,06H15,06L14,08A15,08A19,08J16,08K18,09G15,10G17,11D16,11D17,11H16,11L15,12E15,12K16,13B18,13F16,13F18,13K15,13K17,14C18,15J17,15M15,16F14,16G18,16J15,17E16,17H18,18K18,18L15,20K14,20K16,21A19,21D16,22A16,22F17,22G14,22L18,22M16,23B15,23E16,23K17,25E18,25G17,25J18,26C18,26E15,26H16,26J17,26L15,28C17,28G16,28K15,29B16,29H17,29K18,30J16,31C15, |
| XLT-2220 | 01B17,03B15,03E18,05C15,05K15,06L14,09K17,16L15,19A15,20F14,23B16,23F17,30D15,31C15, |
| XLT-2408 | 03A18,03J14,04G18,07L16,09B15,11G14,14C18,18L14,26A17,27L18,28K14,28K15, |
| XLT-2410 | 09F17,13G16,16B15,24J18,25H15,29F16,29G14, |
| XLT-2412 | 01L18,08J16,09F17,15D16,18H16,18J14,24J18,26C19,27C17,27H15,03G14, |
| XLT-2414 | 01C17,03D17,04E18,06H15,07L18,08A19,09E17,09M15,11B16,11G18,11H16,15K18,16J16,21A19,21D16,23G14,27G18,27L18,29C16,29H17,29K18,30A15,31K17, |

| | |
|----------------------------------|--|
| XLT-2416 | 01B17,01F18,02D19,02E14,02E16,03E17,06C15,06C19,06D18,08A19,09B15,09M15,10B15,10B16,11G14,11L15,12E14,13K15,13M18,14A19,16F17,17E18,18K18,19C15,20C18,20D15,20L18,21A15,22G14,23E17,23G15,23G18,24B16,24L17,25F14,27F18,27L14,28C17,28C19,28K14,29C16,30A18, |
| XLT-2418 | 01D15,01E14,02C17,03M15,07F16,09A19,09B15,09C15,09E17,10B16,12E17,18K18,19C15,19H15,19J14,22F15,22F17,22G14,22G15,22H14,23B16,26C18,26C19,27G18,29C16,29H17,30K17,07H14, |
| XLT-2420 | 03A19,03E16,04G18,07C19,08G15,10L14,12B16,12E17,13K14,17K14,18K18,19H16,19J17,19L14,23G14,25A18,26L14,26L15,29L17,30A18,31C15, |
| X-Ray Stomach Tubes (XST) | |
| XST-2408 | 01D15,01E15,06K15,08M15,09B15,12F17,18B16,18H14,19G17,22K18,25J14, |
| XST-2410 | 19G17,26L15,29A16, |
| XST-2412 | 26L15,27L14, |
| XST-2414 | 01C17,07K15,14K16,29F17, |
| XST-2416 | 04G17,06D17,13K15,18H16, |
| XST-2418 | 02C17,06K15,19G17,22H16, |
| XST-2418/BULK | 01B18,10G17,16H18,19K15,20M16,21G15,25H17,26K18,29G16,29H18, |
| XST-2420 | 05C15, |

Addendum B

Field Safety Notice Acknowledgement Form

| | |
|-------------------|--|
| Name | |
| Position | |
| Company/Institute | |
| Address | |

Confirmation of having none of the products referenced in this FSN.

| | |
|---|--|
| We do not have any stock of the affected product. (✓ tick if this is your answer) | |
|---|--|

| Product Code and Quantity of Affected Product in Stock | | |
|--|---------------|------|
| Product Code | LOT Number(s) | Qty. |
| | | |
| | | |
| | | |
| | | |
| | | |

NB: Please contact your local distributor representative or Pennine directly to confirm quantity and type of quarantined stock before completing returns or plans for destruction. **Evidence will be required prior to and after destruction.**

| | |
|--|--|
| I hereby confirm that I have read and understood the contents of this Field Safety Notice. | |
| Sign | |
| Date | |

Please return the completed reply form to: recalls@penninehealthcare.co.uk